

ARIANA Patient Transfer System:

The board with sliding sleeve for lateral transfers Ariana-Roller can comfortably transfer patients between beds, stretchers, operating tables, X-ray tables, etc. Patients may transfer up to 135 kg between two surfaces with a gap of 20 cm and spaced 20 cm. With improved core and replaceable outer sheath with very low friction and easy sanitation.

- New Core
- X-ray Transparent
- Less Weight
- More Comfortable
- Size:
- Full Board: 180*50 cm
- Half Board: 80*40 cm

R.code	Size
ARM 4011	40*80
ARM 4010	50*180



رول جابجایی بیمار مدل ARIANA:

تخته انتقال بیمار مدل آریانا نسل جدید سیستم های ترانسفر برای جابجایی بیمار به پهلو (در حالت Lateral) از روی تخت بستری، برانکارد، تخت جراحی و یا تخت رادیولوژی (X-ray) می باشد. این تخته قابلیت جابجا کردن بیمار تا ۱۳۵ کیلوگرم بین دو سطح با فاصله حداکثر ۲۰ سانتی متر را دارا می باشد. مواد ارتقا یافته بدنه بسیار سبک بوده و رویه خارجی قابل تعویض آن با اصطکاک بسیار پایین از مزایای قابل توجه این رول می باشد. همچنین رویه رول به راحتی تمیز و در صورت نیاز ضد عفونی می شود. این رول قابل استفاده در بخش رادیولوژی بیمارستان ها بوده و نسبت به X-ray شفاف است و هیچ ردی از خود به جا نمی گذارد.



Step 1



Step 2



Step 3

