


Dual-Sided Gel-Memory Foam Wheelchair Positioning Cushion

High-reticulated viscoelastic gel : avoids gel migration thus preventing all punching of the support. the addition of a high resilience foam base enables the cushion to adept better to bony protuberances, thus providing a larger contact surface, therefore more effective distribution of pressures. Moreover, it improves patient comfort and stability. Because the two layers of the material are moulded together, they are completely bonded, thus eliminating all risks of separating over time.

PRODUCT CODE

AWGC202



Gel (2.5Kg)


ژل

Memory Foam

معموری فوم

Gel (2.5Kg)

ژل



Dual-Sided Gel-Cold Foam Wheelchair Positioning Cushion


High-reticulated viscoelastic gel avoids gel migration thus preventing the punching of the support. The viscoelastic fluid gel ensures reduction of the the Max Pressure on sensitive areas.

The addition of a high-resilience foam base enables the cushion to adept better to bony protuberances, thus providing a larger contact surface therefore more effective distribution of pressure. Moreover it improves patient comfort and stability.

Because the two layers of the material are moulded together, they are completely bonded, thus eliminating all risks of separating over time.

PRODUCT CODE

AWGC201



Gel (1 Kg)


ژل

Cold Foam

فوم سرد

Gel (1 Kg)

ژل



Single-Sided Gel-Cold Foam Wheelchair Positioning Cushion

The Ariana Gel-Foam Wheelchair Cushion consists of a layer of gel and a layer of cold foam. The addition of a high-resistance foam base enhances the cushion's ability to conform to bony prominences, thereby providing a larger contact surface. This leads to more effective pressure distribution

PRODUCT CODE

AWGC101

